

NAVOFIL®

A L W A Y S F O R Y O U R S M I L E

Etchant Gel

Scope

In intraoral prosthesis or orthodontic treatment, etchant gel is used to treat the surfaces of teeth, metal, ceramic and other prostheses to remove the contaminated layer, rough surface and improve their surface properties.

Specifications

2.5ml/syringe

Model

Model	H ₃ PO ₄
AT-20	20%
AT-33	33%
AT-35	35%
AT-37	37%

CE 0197

ATBIO
AT&M Biomaterials Co., Ltd.

☎ 010-56764321

🌐 www.atmbio.com

✉ info@atmbio.com

📍 No.12 Yongcheng Bei Road, Zhongguancun Yongfeng Industry Base, Haidian District, 100094 Beijing, P.R. China.



NANOFIL[®]

A L W A Y S F O R Y O U R S M I L E

GEL-LIKE CONSISTENCY

Preventing unintended flows to reduce risk

EASY TO OPERATE

Easy to be coated

EASY TO CLEAN

Hydrophilic formula for residue-free cleaning

DARK BLUE DESIGN

Easy to identify, no dead space for cleaning

SHORT TIME FOR ETCHING

Rapid surface treatment using phosphoric acid as an effective acid etching factor.

VARIOUS CONCENTRATIONS

A wide range of concentrations to suit different habits or needs



CE 0197

ATBIO

AT&M Biomaterials Co., Ltd.



010-56764321



www.atmbio.com



info@atmbio.com



No.12 Yongcheng Bei Road, Zhongguancun Yongfeng Industry Base, Haidian District, 100094 Beijing, P.R. China.