

NANOFIL[®]

A L W A Y S F O R Y O U R S M I L E

DC Core Composite

Scope

- Permanent cementation for all types of root canal posts
- Core build-ups
- Permanent cementation of crowns, bridges, inlays, onlays (ceramic, metal, and composite)

Specifications

8g/syringe (With 5 mixing heads)

Color

Pale yellow



☎ 010-56764321

🌐 www.atmbio.com

✉ info@atmbio.com

📍 No.12 Yongcheng Bei Road, Zhongguancun Yongfeng Industry Base, Haidian District, 100094 Beijing, P.R. China.

CE 0197

ATBIO

AT&M Biomaterials Co., Ltd.

NAVOFIL[®]

A L W A Y S F O R Y O U R S M I L E

DOUBLE CURING

Self curing and light curing for more complete curing

IT'S VERY ADAPTABLE

All crown and core reconstructions and post stud bonding

GOOD MECHANICAL PROPERTIES

Flexural strength ≥ 130 mpa,
compression strength ≥ 260 mpa

TWO-COMPONENT DESIGN


The intraoral injection head can be operated more precisely in the mouth

FAST CURING SPEED

It will completely solidify after 20 seconds of light

STRONG ADHESIVE PROPERTY

It can be used with quartz pile, metal pile, etc.



010-56764321



www.atmbio.com



info@atmbio.com



No.12 Yongcheng Bei Road, Zhongguancun Yongfeng Industry Base, Haidian District, 100094 Beijing, P.R. China.



ATBIO

AT&M Biomaterials Co., Ltd.